



PRODUCT SPECIFICATION

catalog number:	J65200
product description:	Sex Hormone Binding Globulin (SHBG) Antigen
lot number:	[lot specific]
source:	Normal human serum/plasma
form:	White lyophilized solid; 0.2µm filtered
purity:	>90%
buffer:	Lyophilized from 0.02M NH ₄ HCO ₃ . Buffer salts may be present.
preservative:	none
concentration:	[lot specific] mg/mL
bioburden:	< 100 CFU/mL
expiration date:	[lot specific]
recontrol/extended expiration date:	[lot specific]
reconstitution:	Recommend reconstitution using a phosphate buffer, pH >7.0 containing 0.15M NaCl.
storage:	2-8°C (shipping); upon arrival store at <-60 long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard associated with this product. Donor material has been tested and found negative for HBsAg, anti-HIV 1&2 and anti-HCV by FDA approved methods. However, this material should still be handled with caution as potentially hazardous.