



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J66030**

PRODUCT DESCRIPTION:	Human Parathyroid Hormone (PTH) 1-84, Recombinant Human Antigen
SOURCE:	<i>Escherichia Coli</i>
PURITY:	>95% by SDS-PAGE
FORM:	Lyophilized powder. Each 1.3 mg of lyophilized protein contains 50 µg mannitol and 20 mM PBS, pH 7.0
RECONSTITUTION:	Reconstitute in sterile deionized water, add 0.1% BSA or HSA to improve stability.
STORAGE TEMPERATURE:	Stable at room temperature for 3 weeks, but should be stored at -20°C. Upon reconstitution PTH should be stored at 2-8°C for 2-7 days and long term at -20°C. Prevent freeze-thaw cycles.
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.