



PRODUCT SPECIFICATION

CATALOG NUMBER: J66130

PRODUCT DESCRIPTION:	Ca 19-9 GI Antigen, Cell Culture, High Pure
SOURCE:	Human adenocarcinoma cell line supernatant
BUFFER:	PBS Buffer, pH 7.4±0.2, 0.2 µ filter
PRESERVATIVE:	0.1% sodium azide
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
CROSS REACTIVITY:	Ca 15-3 ≤ 25% Ca 125 ≤ 25%
EXPIRATION DATE:	Three years from date of manufacture
BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. The cell line supernatant was tested and found negative for anti-HIV-1/HIV-2, HBsAg, anti-HCV, HIV-1 and HCV by Nucleic Acid Test per current FDA methods.