



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J66202**

PRODUCT DESCRIPTION:	Ca 125 Ovarian Cancer Antigen, Cell Culture, Part Pure
SOURCE:	Cell culture supernatant
PURITY:	≥1600 Ku activity/ml per OD at 280 nm
BUFFER:	PBS Buffer, 0.2 μ filter
PRESERVATIVE:	0.1% sodium azide
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
CROSS REACTIVITY:	Ca 15-3 ≤ 25% Ca 19-9 ≤ 25%
EXPIRATION DATE:	Two and a half years from date of manufacture
BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. The cell line supernatant was tested and found negative for anti-HIV-1/HIV-2, HIVag, HBsAg, and HCV.