



PRODUCT SPECIFICATION

CATALOG NUMBER: J66210

PRODUCT DESCRIPTION:	CA 125, Ovarian Cancer Antigen
SOURCE:	Human fluid
BUFFER:	25mM Tris, 140 mM NaCl, pH 7.4
PRESERVATIVE:	0.15% sodium azide
FINAL FILTRATION:	0.2 Micron filtered
RELATIVE PROTEIN:	36714.66 units/mL/OD
TOTAL PROTEIN:	5.765 OD units at 280 nm
CROSS REACTIVITY:	CA 19-9 6.47% CA 15-3 1.39% CEA 37.87 ng/mL AFP 0.0 Transferrin 0.0
STORAGE TEMPERATURE:	2-8°C short term
BIOHAZARD:	All donor units have been tested and found negative for HIV 1/2 and HCV antibodies and HBsAg. This material contains sodium azide as a preservative, take appropriate precautions.