



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J66301**

PRODUCT DESCRIPTION:	Ca 15-3 Breast Tumor Antigen, Cell Culture, High Pure
SOURCE:	Human BTA cell line supernatant
PURITY:	>5 Ku activity/ml per OD at 280 nm
BUFFER:	PBS, pH 7.4±0.2, 0.2 µ filter
PRESERVATIVE:	0.1% sodium azide
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
CROSS REACTIVITY:	CA 125 ≤ 25% CA 19-9 ≤ 25%
EXPIRATION DATE:	Two years from date of manufacture
BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. The cell line supernatant was tested and found negative for HbsAg, anti-HCV and anti-HIV 1 & 2 antibodies; HIV-1 and HCV by Nucleic Acid Test per current FDA methods.