



PRODUCT SPECIFICATION

catalog number:	J66330
product description:	CA 15-3 Breast Cancer Antigen, Standard Grade
lot number:	[lot specific]
source:	Human fluids
form:	Liquid
buffer:	1% Tritron X100, 1M NaCl
preservative:	0.1% NaN ₃
concentration:	[lot specific] KU/mL
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	2-8°C short term (up to 12 months); -20°C long term (up to 36 months)
hazard/biohazard:	<p>This material contains sodium azide as a preservative; take appropriate precaution.</p> <p>Donor material has been tested and found negative for HBsAg, anti-HCV and anti-HIV 1&2 by FDA approved methods. However, this material should still be handled as potentially hazardous.</p>