



## PRODUCT SPECIFICATION

catalog number:	<b>J75000</b>
product description:	<i>H. pylori</i> Antigen
lot number:	[lot specific]
source:	Derived from one single <i>H. pylori</i> colony; Cag+.
form:	Detergent extract of whole <i>H. pylori</i> cell sonicates.
buffer:	0.015M KPO <sub>4</sub> , 0.85% NaCl, pH 7.4
preservative:	0.1% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL (by Pierce BCA protein assay)
specificity:	<i>H. pylori</i> specific
storage:	-20°C or below
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution.  This material is extracted from whole <i>H. pylori</i> cells and should be handled with caution as potentially hazardous.