



## PRODUCT SPECIFICATION

catalog number:	<b>J75110</b>
product description:	<i>H. pylori</i> Antigen, Part Pure
lot number:	[lot specific]
source:	<i>In vitro</i> purified proteins from <i>H. pylori</i> ; <i>Urease+</i> .
form:	liquid
buffer:	15mM phosphate, 0.15M NaCl, pH 7.0
preservative:	0.1% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.