



PRODUCT SPECIFICATION

CATALOG NUMBER: J80110

PRODUCT DESCRIPTION:	NSE Recombinant Antigen, High Pure
SOURCE:	<i>E. coli</i>
FORM:	Liquid
APPEARANCE:	Clear to colorless
PURIFICATION:	Metal affinity chromatography
PURITY:	>95% by SDS-PAGE ^{12%} Major band @ 50.8 kD
BUFFER:	Aseptically filtered liquid in 20mM sodium phosphate, 150mM NaCl, 50% Glycerol, pH 7.4 ± 0.2
BIO-BURDEN:	< 10 CFU/ml
STORAGE TEMPERATURE:	-20°C. Avoid multiple freeze thaw cycles.
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.