



PRODUCT SPECIFICATION

CATALOG NUMBER: J80250

PRODUCT DESCRIPTION:	Human Serum Albumin Antigen, Low Folate
SOURCE:	Human plasma
FORM:	Powder
PURITY:	>98%
STORAGE TEMPERATURE:	2-8°C
ALBUMIN:	99.3%
PH:	6.9
SODIUM:	30.6 mEq/L
CHLORIDE:	8 mEq/L
POTASSIUM	0.07 mEq/L
CALCIUM:	3.5 mg/dL
OPTICAL DENSITY:	<u>0.05@540nm</u>
IRON:	64 ug/dl
AMMONIA:	9 umol/L
MICROBIAL(7%)	0.0 cfu/ml
VITAMIN B12:	158 pg/mL
FOLATE:	1.4 ng/mL
MOISTURE:	1.8%
BIOHAZARD:	Donor plasma has been tested and found non reactive for HIV 1&2 antibody, HCV antibody, and HBsAg.