



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: J80400**

PRODUCT DESCRIPTION:	Apolipoprotein A1
SOURCE:	Human plasma
FORM:	Liquid frozen in 0.1M NH <sub>4</sub> HCO <sub>3</sub>
PURITY:	Nominally greater than 95% by SDS PAGE
STORAGE TEMPERATURE:	-15°C
BIOHAZARD:	Starting material donor tested and found negative for HBsAg and antibodies to HIV 1&2 and HCV. There is no known biohazard associated with this product.