



PRODUCT SPECIFICATION

CATALOG NUMBER: J81000

PRODUCT DESCRIPTION:	Human Complement C3
SOURCE:	Normal human serum
FORM:	Single homogenous batch supplied lyophilized
PURITY:	96%
BUFFER:	sodium phosphate buffer, pH 7.2, containing 0.15M NaCl
PRESERVATIVE:	No preservative added
STORAGE TEMPERATURE:	2-8°C
BIOHAZARD:	This product has been tested at serum donor level by FDA approved methods and found negative for HIV I & II antibodies and Hepatitis B surface antigen and Hepatitis C antibodies.