



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J81810**

PRODUCT DESCRIPTION:	Ferritin
PURITY:	>96%
BUFFER:	0.05M Tris buffer, pH 7.5, 1.0M NaCl
PRESERVATIVES:	0.09% NaN <sub>3</sub>
STORAGE TEMPERATURE:	2-8°C
TESTING:	Donor tested negative HIV 1&2 antibodies, hepatitis B surface antigen, Hepatitis C antibodies, and HIV antigen by FDA approved methods.
ANALYSIS:	Polyacrylamide gradient gel electrophoresis shows one major band.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.