



PRODUCT SPECIFICATION

CATALOG NUMBER: J81840

PRODUCT DESCRIPTION:	Human Placental Ferritin, High Pure
SOURCE:	Human placenta
PURITY:	> 96% by SDS-PAGE ^{7%}
BUFFER:	50 mM Tris-HCl, 150 mM NaCl, pH 7.4
PRESERVATIVE:	0.09% NaN ₃
BIO-BURDEN:	< 10 CFU/ml
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	Donor starting material tested and found negative for HBsAg, HIV 1&2 antibodies and HCV antibody by FDA approved methods. In addition, tested negative for HBV, HIV I, and HCV viral DNA/RNA by PCR. However, this material should still be handled as potentially hazardous. This material contains sodium azide as a preservative, take appropriate precaution.