



PRODUCT SPECIFICATION

catalog number:	J81850428
product description:	Ferritin Human Recombinant, Light Chain
lot number:	[lot specific]
source:	E. Coli
molecular weight:	20 kDa (single, non-glycosylated polypeptide chain 175 amino acids)
form:	Liquid, sterile filtered
purity:	≥ 90% by SDS-PAGE
buffer:	20 mM Tris-HCl, pH 7.5
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C for 2-4 weeks; -20°C long term; add 0.1% HSA or BSA as carrier protein for long term storage; avoid multiple freeze/thaw cycles
hazard/biohazard:	No known biohazards exist for this product