



#### PRODUCT SPECIFICATION

catalog number:	<b>J81851428</b>
product description:	Ferritin Human Recombinant, Heavy Chain
lot number:	[lot specific]
source:	E. Coli
molecular weight:	21 kDa (single, non-glycosylated polypeptide chain 183 amino acids)
form:	Liquid, sterile filtered
purity:	≥ 95% by SDS-PAGE
buffer:	20 mM Tris-HCl, pH 7.5, 1mM DTT with 10% glycerol
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C for 2-4 weeks; -20°C long term; add 0.1% HSA or BSA as carrier protein for long term storage; avoid multiple freeze/thaw cycles
hazard/biohazard:	No known biohazards exist for this product