



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: J82100**

PRODUCT DESCRIPTION:	Human IgA Antigen
SOURCE:	Pooled normal human serum
FORM:	Sterile-filtered liquid
BUFFER:	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
PRESERVATIVE:	None
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.
COUNRTY OF ORIGIN	USA