



PRODUCT SPECIFICATION

catalog number:	J82320
product description:	IgE Antigen
lot number:	[lot specific]
source:	Cell culture supernatant
form:	Liquid, 0.2µm filtered
purity:	[lot specific] %, part pure
buffer:	0.1M Tris, 0.2M NaCl, pH 7.5±0.1
preservative:	none
concentration:	[lot specific] mg/mL (cobas IGE assay)
activity:	[lot specific] KU/mL (Centaur IGE assay)
bioburden:	0 CFU/mL
total protein:	[lot specific] OD units @ A ₂₈₀
manufacture date:	[lot specific]
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. Donor material tested and found non-reactive for HBsAg, anti-HCV, anti-HIV 1&2 and HIV-1 and HCV by NAT per FDA approved methods. However, this material should still be handled with caution as potentially hazardous.