



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J82400**

PRODUCT DESCRIPTION:	Purified Human, IgG antigen
FORM:	Liquid
PURITY:	> 95% by SDS-PAGE
BUFFER:	10 mM Sodium Phosphate, 0.85% (w/v) Sodium Chloride, pH 7.2 (PBS)
PRESERVATIVES:	0.05% NaN <sub>3</sub>
STORAGE TEMPERATURE:	2-8°C short term, -20°C or colder long term
TESTING:	Donor plasma tested and found negative for HBsAg, HIV-1, HIV-2, and HCV antibodies by an FDA approved tests.