



PRODUCT SPECIFICATION

CATALOG NUMBER: J82410

PRODUCT DESCRIPTION:	Human IgG, Fc fragment
SOURCE:	Normal human Serum
FORM:	Liquid
PURIFICATION:	Fractionation and DEAE chromatography. Fc fragments were produced by papain digestion and affinity chromatography.
BUFFER:	Phosphate Buffered Saline, pH 7.2
PRESERVATIVE:	0.09% NaN ₃
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.
HAZARD TESTING:	Donors have been tested and found negative for HIV-1 antigen and antibodies for HIV and HCV and non-reactive for HBsAg.