



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J82420493**

PRODUCT DESCRIPTION:	Human IgG F(ab)'2 fragment
LOT NUMBER:	Reference
FORM:	Liquid, 0.2 µm filtered
PURIFICATION:	Protein A chromatography
BUFFER:	10 mM sodium phosphate, 0.15 sodium chloride, pH 7.2
PRESERVATIVE:	0.05% (w/v) sodium azide
CONCENTRATION:	8.32 mg/ml (E <sup>1%</sup> at 280 nm = 13.0)
STORAGE TEMPERATURE:	2-8°C short term
HAZARD/BIOHAZARD:	Donor serum was tested and found negative for antibodies to HIV 1/2, HCV, HBsAg, HIV-1 RNA, HCV RNA and Syphilis by FDA approved methods. This material contains sodium azide as a preservative, take appropriate precaution.