



## PRODUCT SPECIFICATION

catalog number:	<b>J82820</b>
product description:	Lambda Free light Chain, High Pure
lot number:	[lot specific]
source:	Pooled human urine
form:	Liquid, sterile filtered
buffer:	50 mM Sodium Phosphate, 0.15M NaCl, pH 7.4 ± 0.2
preservative:	0.1% Sodium Azide
purity:	>95% by SDS PAGE (reducing conditions)
protein activity:	[lot specific] mg/mL (Binding Site Nephelometric method)
cross reactivity:	Kappa Free light chain ≤1% IgG ≤1% IgA ≤1% IgM ≤1% Cystatin C ≤1% Microalbumin ≤1%
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	Negative for HBsAg, Anti-HCV and Anti-HIV I & II , syphilis by CE Marked test method