



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J83800**

PRODUCT DESCRIPTION:	Human Thyroid Peroxidase Antigen
SOURCE:	Human thyroid gland
FORM:	Lyophilized
PURITY:	>90% pure by SDS-PAGE
RECONSTITUTION:	Use of distilled water is recommended.
STORAGE:	2-8°C for shipping purposes, -20°C long term in lyophilized form and after reconstitution
HAZARD/BIOHAZARD:	Donors have been tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV antibody, and syphilis.