



PRODUCT SPECIFICATION

catalog number:	4-001-A (J84400)
product description:	Purified Goat IgG
lot number:	[lot specific]
source:	Normal goat serum
form:	Liquid
purity:	> 95% pure by SDS-PAGE
buffer:	10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% (w/v) sodium azide
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C
country of origin:	USA
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.