



## PRODUCT SPECIFICATION

catalog number:	<b>J91009</b>
product description:	Human Procalcitonin (PCT) Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i>
purification:	Ion exchange chromatography
form:	Lyophilized
purity:	≥97% by SDS-PAGE 15% (Major band @ 13 kDa)
buffer:	Lyophilized from 20mM Sodium Phosphate, 150mM NaCl, 0.1M Lysine, 5% Trehalose, pH 7.4 ±0.2
preservative:	Proclin 300
bio-burden:	≤ 10 CFU/ml
storage:	-20°C
hazard/biohazard:	There is no known hazard or biohazard associated with this product.