



PRODUCT SPECIFICATION

catalog number:	J91010
product description:	Human Procalcitonin (PCT) Antigen, Recombinant
lot number:	[lot specific]
source:	<i>E. coli</i>
form:	Liquid
purity:	> 95%
buffer:	20 mM Tris HCl, 50% glycerol, pH8.0
preservative:	none
molecular weight:	~17 kDa
concentration:	[lot specific] mg/mL
storage:	Shipped at 2-8°C; upon arrival, store at -20°C avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.