



PRODUCT SPECIFICATION

| | |
|----------------------|--|
| catalog number: | J91011 |
| product description: | Human Procalcitonin (PCT) Antigen, Recombinant |
| lot number: | [lot specific] |
| source: | Mammalian |
| form: | Liquid |
| purity: | >95% |
| buffer: | 1x PBS, 50% glycerol, pH 7.4 |
| preservative: | none |
| molecular weight: | 13.45 kDa |
| concentration: | [lot specific] mg/mL |
| storage: | Shipped at 2-8°C; upon arrival, store at -20°C avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | There is no known hazard or biohazard associated with this product. |