



PRODUCT SPECIFICATION

CATALOG NUMBER: J91200

PRODUCT DESCRIPTION:	GC Binding Protein (Vitamin D Binding Protein)
SOURCE:	Human plasma
FORM:	Lyophilized powder
BUFFER:	Lyophilized from a 20mM sodium phosphate, pH 7.4 with 150mM NaCl
PURITY:	>95% by SDS-PAGE
STORAGE TEMPERATURE:	-20°C or colder
RECONSTITUTION:	Reconstitute in sterile de-ionized H ₂ O (note: to insure complete recovery of product centrifuge briefly prior to opening).
HAZARD/BIOHAZARD:	Source plasma is shown to be non-reactive for HBsAg, anti-HCV, anti-HBc and negative for anti-HIV 1 & 2 by FDA-approved tests.