



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: L82300**

PRODUCT DESCRIPTION:	IgE (Immunoglobulin E) Myeloma Positive Plasma
LOT NUMBER:	Reference
IgE CONTENT:	6,180,000 IU/MI
PRESERVATIVE:	No preservative added
FORM:	Liquid
STORAGE TEMPERATURE:	-20°C or below
HAZARD/BIOHAZARD:	Each human blood donor unit of the above referenced product has been tested and found negative or nonreactive for HBsAg, HIV 1 / 2 antibody, HCV antibody and STS by an FDA approved method.