



PRODUCT SPECIFICATION

CATALOG NUMBER: P11000

PRODUCT DESCRIPTION:	TSH Free Human Plasma
STARTING MATERIAL:	Delipidated, defibrinated human plasma.
PRESENTATION:	Affinity stripped at least 2X to deplete TSH, 0.2u filtered
ANALYTE ASSAY RESULTS:	TSH <0.01 miu/L as determined by Abbott Architect
BIOBURDEN:	<10 cfu/mL
STORAGE TEMPERATURE:	Store below -15°C; shelf life minimum 1 year for HIV I & II antibodies, Hepatitis B surface antigen, and Hepatitis C antibodies by FDA approved methods.
HAZARD/BIOHAZARD:	Starting material donor tested and found negative. This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.