



PRODUCT SPECIFICATION

CATALOG NUMBER: P19000

PRODUCT DESCRIPTION:	T3, T4 Free Human Serum, Lipid Stripped	
PRESERVATIVES:	None	
TOTAL PROTEIN:	6.3 gm/dl	
STORAGE TEMPERATURE:	-20°C	
BIOCHEMICAL ASSAY:	Albumin	4.3 gm/dl
	Globulin	2.0 gm/dl
	Cholesterol	<3.0 mg/dl
	Triglycerides	10.0 mg/dl
	Uric Acid	<0.20 mg/dl
	Glucose	<2.0 mg/dl
	Creatinine	<0.20 mg/dl
	Calcium	<0.2 mg/dl
	Phosphorus	1.9 mg/dl
	Total Bilirubin	0.10 mg/dl
	Sodium	132.0 mEq/L
	Potassium	1.88 mEq/L
	Chloride	92.1 mEq/L
	ALT	0.0 U/L
	AST	<4.0 U/L
	GGT	11.0 U/L
	Bun	<2.0 mg/dl
	Alkaline Phosphatase	16.0 U/ml
BIOBURDEN:	<100 cfu/mL	
T3/T4 ASSAY RESULTS:	<Sensitivity level of assays used for measurement	
NOTE:	Donors tested and found non-reactive for HBsAg, anti-HIV 1/2, anti-HCV, and HIV-1 Ag by FDA approved methods.	