



PRODUCT SPECIFICATION

CATALOG NUMBER: P21400

PRODUCT DESCRIPTION:	Cortisol Stripped Serum	
PRESERVATIVES:	None	
TOTAL PROTEIN:	6.2 gm/dl	
STORAGE TEMPERATURE:	-20°C long term	
BIOCHEMICAL ASSAY:	Albumin	4.3 gm/dl
	Globulin	1.9 gm/dl
	Cholesterol	2.0 mg/dl
	Triglycerides	<4.0 mg/ml
	Total Bilirubin	<0.1 mg/dl
	Uric Acid	0.0 mg/dl
	Glucose	<2.0 mg/dl
	Creatinine	0.0 mg/dl
	ALT	3.0 U/L
	AST	14.0 U/L
	GGT	13.0 U/L
	Alkaline Phos.	18.0 U/L
	Calcium	<0.2 mg/dl
	Phosphorus	2.3 mg/dl
	Sodium	133.7 mEq/L
	Potassium	2.15 q/L
	Chloride	92.7 mEqL
ANALYTE ASSAY RESULTS:	Cortisol <0.02 ng/ml	
CFU/ML AT BOTTLING:	<100	
HBsAg Antigen:	Tested negative	
HIV1&2, HCV Antibody:	Tested negative	
STS	Tested negative	
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.	