



PRODUCT SPECIFICATION

CATALOG NUMBER: P30000

| | |
|----------------------|---|
| PRODUCT DESCRIPTION: | IgG Free Human Serum, Lipid Stripped |
| PRESERVATIVES: | None |
| TOTAL PROTEIN: | 6.0 gm/dL |
| BIOCHEMICAL ASSAY: | Albumin 4.8 gm/dL Globulin 1.2 gm/dL Cholesterol 1.0 mg/dL Triglycerides <4.0 mg/dL TotalBilirubin <0.1 mg/dL Glucose <2.0 mg/dL Creatinine 0.0 mg/dL ALT 6.0 U/L AST 13.0 U/L GGT 16.0 U/L Alk. Phos. 1.0 U/L Calcium <0.2 mg/dL Inorganic Phosphate 1.6 mg/dL Sodium 145.7 mEq/L Potassium 1.6 mEq/L Chloride 103.4 mEq/L Uric Acid 0.0 mg/dL |
| MICROBIOLOGY ASSAY: | <100 CFU/mL at bottling |
| STORAGE TEMPERATURE: | -20°C |
| VIRAL TESTING: | Donor units are tested and found non-reactive for HIV 1&2 antibody, HCV antibody, HTLV-I/2 antibody, HBsAg, STS, and HIV-1. |