



## PRODUCT SPECIFICATION

**CATALOG NUMBER: P61000**

PRODUCT DESCRIPTION:	AFP Free Human Serum, Lipid Stripped
PRESERVATIVES:	None
TOTAL PROTEIN:	5.3 gm/dl
BIOCHEMICAL ASSAY:	Albumin 3.5 gm/dl Globulin 1.8 gm/dl Cholesterol 1.0 mg/dl Triglycerides <4.0 mg/dl Uric Acid 0.0 mg/dl Glucose <2.0 mg/dl Creatinine 0.0 mg/dl Calcium <0.2 mg/dl Total Bilirubin <0.1 mg/dl Sodium 133.6 mEq/L Potassium 2.17mEq/L Chloride 89.80 mEqL ALT 5.0 U/L AST 2.0 U/L GGT 8.0 U/L Phosphorus 2.1 mg/dl Alkaline Phosphatase 24.0 U/ml
BIOBURDEN:	<100 cfu/mL
T3/T4 ASSAY RESULTS:	<Sensitivity level of assays used for measurement
STORAGE TEMPERATURE:	-20°C
NOTE:	Donors tested and found negative for HBsAg, anti-HIV 1/2, anti-HCV, and HIV-1 Ag by FDA approved methods.