



PRODUCT SPECIFICATION

CATALOG NUMBER: P62000

PRODUCT DESCRIPTION:	CEA Free Serum, Lipid Stripped	
PRESERVATIVES:	None	
FILTRATION:	0.2 micron filtered	
TOTAL PROTEIN:	6.2 gm/dL	
CEA CONTENT:	Below the detectable limit	
CHEMICAL ASSAY:	Glucose	<2.0 mg/dL
	Creatinine	0.0 mg/dL
	Albumin	4.5 gm/dL
	Globulin	1.7 gm/dL
	Total Bilirubin	0.15 mg/dL
	ALT	9.0 U/L
	AST	16.0 U/L
	GGT	17.0 U/L
	Alk. Phos.	34.0 U/L
	Sodium	139.2 mEq/L
	Potassium	1.84 mEq/L
	Chloride	96.1 mEq/L
	Calcium	<0.2 mg/dL
	Phosphorus	1.4 mg/dL
	Uric Acid	0.0 mg/dL
	Cholesterol	2.0 mg/dL
	Triglycerides	<4.0 mg/dL
MICROBIOLOGY ASSAY:	Colony forming units/mL at bottling <100	
STORAGE TEMPERATURE:	-20°C	
NOTE:	Each human blood donor Unit of the above referenced product has been tested by FDA approved methods and found to be non reactive for HIV-1/2 Antibody, HCV antibody, HTLV-1/2 Antibody, Syphilis (STS) and HBsAg. The above referenced product was only formed with these donors.	