



PRODUCT SPECIFICATION

CATALOG NUMBER: P81000

PRODUCT DESCRIPTION: Ferritin Free Human Serum, Lipid Stripped

PRESERVATIVES: None

TOTAL PROTEIN: 6.1 gm/dl

FERRITIN CONTENT: <1.0 ng/ml by Elecsys

BIOCHEMICAL ASSAY:

Albumin	4.1 gm/dL
Globulin	2.0 gm/dL
Cholesterol	1.0 mg/dL
Triglycerides	<4.0 mg/dL
Total Bilirubin	<0.1 mg/dL
Uric Acid	0.0 mg/dL
Glucose	<2.0 mg/dL
Creatinine	0.0 mg/ml
Phosphorus	2.1 mg/dl
Alk Phos	47.0 U/L
Calcium	<0.2 mg/ml
Sodium	139.0 mEq/L
Potassium	2.11 mEq/L
Chloride	97.7 mEq/L
ALT	6.0 U/L
AST	14.0 U/L
GGT	11.0 U/L

MICROBIOLOGY ASSAY: CFU/mL <100 at bottling

STORAGE TEMPERATURE: -20°C

BIOHAZARD: Each human blood donor Unit of the above referenced product has been tested by FDA approved methods and found to be non reactive for HIV-1/2 Antibody, HCV antibody, HIV Ag, and HBsAg. The above referenced product was only formed with these donors.