



## PRODUCT SPECIFICATION

**CATALOG NUMBER: T82422**

PRODUCT DESCRIPTION:	Affinity Purified Goat Anti Human IgG-Fc Alkaline Phosphatase Conjugate
FORM:	Liquid
BUFFER:	1.0MM MgCl <sub>2</sub> , 0.1MM ZnCl <sub>2</sub>
STORAGE TEMPERATURE:	2-8°C
COUNTRY OF ORIGIN:	USA
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.