



## PRODUCT SPECIFICATION

**CATALOG NUMBER: T82451**

PRODUCT DESCRIPTION:	Affinity Pure Goat anti Human IgG H&L, FITC Conjugate (cross absorbed using bovine, chicken, goat, horse, pig, rabbit, and rat immunosorbents)
FORM:	Liquid
BUFFER:	PBS, pH 7.2
PRESERVATIVES:	0.1% Sodium Azide
STABILIZER:	0.2% BSA
STORAGE TEMPERATURE:	2-8°C