



PRODUCT SPECIFICATION

CATALOG NUMBER: T82730

PRODUCT DESCRIPTION:	Affinity Pure Goat anti Human IgA, F(ab') ₂ HRP conjugate
FORM:	Freeze dried powder
BUFFER:	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
PRESERVATIVE:	0.01% thimerosal
SUGGESTED DILUTION:	1/500-1/5,000 for immunohistochemistry on tissue secretions, 1/5,000-1/100,000 for ELISA and Western blots
STORAGE:	4-8°C short term, -20°C long term