



PRODUCT SPECIFICATION

catalog number:	V19030
product description:	Triiodothyronine (T3)-HRP Conjugate
lot number:	[lot specific]
buffer:	0.02M phosphate, pH 7.5, 1:1 glycerol
concentration:	[lot specific] mg/mL; HRP 250 U/mg
titer:	1:10,000
position:	T3
storage:	2-8°C short term (shipping); upon arrival store immediately at -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.