



PRODUCT SPECIFICATION

| | |
|----------------------|----------------------------------------------------------------------------------------------------------------------|
| catalog number: | V19110 |
| product description: | Thyroxine (T4)-ALP Conjugate |
| lot number: | [lot specific] |
| buffer: | 0.005M Tris HCl, 0.005M MgCl ₂ , 0.0001M ZnCl ₂ , pH=7 50% Glycerol |
| concentration: | [lot specific] mg/mL |
| storage: | 2-8°C short term (shipping); upon arrival store immediately at -20°C long term avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | There is no known biohazard associated with this product. |