



PRODUCT SPECIFICATION

catalog number:	V50045
product description:	K2-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
buffer:	PBS, pH 7.1
preservative:	none
concentration:	[lot specific] mg/mL
pairing:	This conjugate can be paired with K2 polyclonal antibody D50042, D50043 and D50045.
notes/applications:	This conjugate is known to work well in lateral flow.
storage:	2-8°C (\leq 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.