



PRODUCT SPECIFICATION

catalog number:	V51710
product description:	Zolpidem-HRP Conjugate
lot number:	[lot specific]
form:	Liquid; enzyme stabilizing diluent
buffer:	PBS, pH 7.0
preservative:	0.005% thimerosal and other additives
concentration:	[lot specific] mg/mL
pairing:	This conjugate can be paired for with zolpidem polyclonal antibodies D51710, D51720 or D51730.
notes/applications:	This conjugate is known to work well for ELISA.
storage:	2-8°C
hazard/biohazard:	This material contains thimerosal as a preservative; take appropriate precaution. There is no known biohazard associated with this product.