



## PRODUCT SPECIFICATION

catalog number:	<b>V51720</b>
product description:	Zolpidem-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
buffer:	Distilled water
preservative:	none
concentration:	[lot specific] mg/mL
pairing:	This conjugate can be paired for with zolpidem polyclonal antibodies D51710, D51720 or D51730.
notes/applications:	This conjugate is known to work well for lateral flow assay.
storage:	2-8°C short term ( $\leq$ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.