



PRODUCT SPECIFICATION

catalog number:	V53044
product description:	Fentanyl-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
purification method:	Chromatography purified
purity:	>95%
buffer:	10 mM PBS, pH 7.4
preservative:	0.1% NaN ₃
concentration:	[lot specific] mg/mL
hapten/protein:	~100:1
link position:	Position 4
pairing:	When used as the coating antibody, this conjugate forms a match pair with A530445 (as the label).
notes/applications:	ELISA and Lateral flow
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.