



PRODUCT SPECIFICATION

catalog number:	V56045
product description:	Ethyl Glucuronide (EtG)-BSA Conjugate
lot number:	[lot specific]
form:	Liquid
buffer:	0.01 M sodium phosphate, 0.15M NaCl, pH 7.4
preservative:	none
concentration:	[lot specific] mg/mL
hapten/protein:	16-18
pairing:	This conjugate can be paired with EtG monoclonal antibody A56044.
notes/applications:	This conjugate is known to work well in lateral flow and ELISA.
storage:	2-8°C short term (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.