



## PRODUCT SPECIFICATION

catalog number:	<b>V56048</b>
product description:	Ethyl Glucuronide (EtG)-HRP Conjugate
lot number:	[lot specific]
form:	Liquid; enzyme stabilizing diluent
buffer:	Tris-HCl, pH 7.0
preservative:	none
concentration:	[lot specific] mg/mL
pairing:	This conjugate can be paired with EtG monoclonal antibody A56044.
notes/applications:	This conjugate is known to work well in ELISA.
storage:	2-8°C short term ( $\leq$ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	There is no known hazard or biohazard associated with this product.