



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>V91142</b>   |
| product description: | <b>Vitamin D Derivative – Carboxylic Acid</b>                       |
| lot number:          | [lot specific]  |
| form:                | Lyophilized   |
| concentration:       | [lot specific] mg/mL  |
| reconstitution:      | Reconstitute with Ethanol   |
| storage:             | -20°C before and after reconstitution in inert gas environment      |
| pairing:             | This antigen can be used with A91156, A91157 or A91158 as conjugate |
| hazard/biohazard:    | There is no known hazard or biohazard associated with this product. |